

Model 1470

Channel Firmware History

Purpose

This document summarizes the changes made to the Channel firmware fitted to the Model 1470E since version 3.12

Channel Firmware History

The table below lists the major firmware revisions. Revision numbers missing from the sequence were not formally released.

Revision	Date	Changes Included
3.13	3 June 2003	Add a new feature, measure-on-change.
3.14	19 September 2003	Fix step duration and the auxiliary scaling bugs.
3.15	30 April 2004	Correct errors in the time stamp.
3.16	25 November 2004	New release for 1470E
3.17	22 December 2009	Implemented Fast/Normal Mode change requested by Scribner
3.18	4 April 2007	Corrects a bug in the interface to the Multistat software from Scribner Associates.
3.19	16 June 2009	Updated to fix a problem where spikes were seen in the acquired data plots when using a booster and changing steps.



...part of **AMETEK**® Advanced Measurement Technology

UNIT B1 ARMSTRONG MALL
SOUTHWOOD BUSINESS PARK
FARNBOROUGH, GU14 0NR
UNITED KINGDOM
Phone: +44 (0) 1252 556 800
Fax: +44 (0) 1252 556 899

801 SOUTH ILLINOIS AVENUE
OAK RIDGE
TN 37831-2011
USA
Phone: +1 865 425 1360
Fax: +1 865 481 2410

Visit our website for a complete list of our global offices and authorized agents

solartron.info@ametek.com

www.solartronanalytical.com