

Model 1470

Temperature Firmware History

Purpose

This document summarizes the changes made to the Temperature firmware fitted to the Model 1470E since version 1.0

Channel Firmware History

The table below lists the major firmware revisions. Revision numbers missing from the sequence were not formally released.

Revision	Date	Changes Included
1.1	20 September 2002	The channel firmware has a new feature, measure-on-change. This requires an initial temperature reading.
1.2	4 March 2009	Updated to fix a problem where the temperature reading gave occasional readings in error by an amount equal to the ambient temperature. The problem happened when the temperature was exactly on a few specific temperature points.



...part of **AMETEK**® Advanced Measurement Technology

UNIT B1 ARMSTRONG MALL
SOUTHWOOD BUSINESS PARK
FARNBOROUGH, GU14 0NR
UNITED KINGDOM
Phone: +44 (0) 1252 556 800
Fax: +44 (0) 1252 556 899

801 SOUTH ILLINOIS AVENUE
OAK RIDGE
TN 37831-2011
USA
Phone: +1 865 425 1360
Fax: +1 865 481 2410

Visit our website for a complete list of our global offices and authorized agents

solartron.info@ametek.com

www.solartronanalytical.com